

# PowerFect™ In Vitro DNA and siRNA Transfection Kit (Ver. II)

----- A General Protocol for Transfecting Mammalian Cell

- 100 µl
- 500 µl
- 1000 µl



9601 Medical Center Drive  
Rockville MD 20850  
FAX. 301-560-4919  
TEL. 301-330-5966  
Toll Free. 1-(866)-918-6812  
Email: [info@signagen.com](mailto:info@signagen.com)  
Web: [www.signagen.com](http://www.signagen.com)

This product is for laboratory research ONLY and not for diagnostic use

## Introduction:

Based on our innovative and proprietary lipid-conjugation technology, PowerFect™ Transfection Kit is a liposome based DNA & siRNA delivery tool which was formulated with our proprietary pH Dependent Conformational Change (PDCC) technology to give efficient and reproducible gene knockdown on variety of mammalian cells. PowerFect™ Transfection Kit is the most powerful yet very gentle gene delivery tool for a variety of applications including plasmid DNA and/or siRNA for most of mammalian cell types. Compared with leading products in the market, PowerFect™ is more cost-effective and always provides higher transfection efficiency with less cytotoxicity.

## Contents Per Kit:

- 1x 1.0 mL of PowerFect™ In Vitro Transfection Reagent
- 1x 8.0 mL of PowerFect™ Transfection Buffer (5x )

## Important Guidelines for Transfection:

- PowerFect™ reagent was formulated for DNA and siRNA transfection. The following standard protocol is given for DNA and siRNA transfection to mammalian cells. For a protocol of siRNA/DNA co-transfection, please email us at [info@signagen.com](mailto:info@signagen.com)
- For better efficiency, choosing PowerFect™ Transfection Buffer working solution (1x ) is a must.
- To lower cytotoxicity, transfect cells in presence of serum (10%) and antibiotics.

## Part I. A General Protocol for DNA Transfection.

### Step I. Preparation of Working Solution of PowerFect™ Transfection Buffer (1x )

PowerFect™ Transfection Buffer (5x ) is provided as 5x concentrated stock solution. To make working solution, dilute one part of the stock solution with 4 parts of ddH<sub>2</sub>O. The PowerFect™ Transfection Buffer (1x ) working solution is stable at RT for 24 months.

**Note:** Always keep PowerFect™ Transfection Buffer (5x ) at RT. If refrigerated, white precipitates may appear. It won't affect the transfection efficiency. After dilution with 4 parts of ddH<sub>2</sub>O to make PowerFect™ Transfection Buffer (1x ) working solution, the white precipitates will disappear. Always keep PowerFect™ Transfection Buffer working solution (1x ) at RT.

### Step II. Cell Seeding:

Cells should be plated 18 to 24 hours prior to transfection so that the monolayer cell density reaches to the optimal ~70% confluency at the time of transfection. Complete culture medium with serum and antibiotics is freshly added to each well 30-60 min before transfection.

**Note:** High serum levels (>5%) with antibiotics do NOT have inhibitory effect on transfection efficiency. We recommend using complete serum/antibiotics-containing medium to grow the cells during transfection.

### Step III. Preparation of PowerFect™ -DNA Complex and Transfection Procedures:

For different cell types, the optimal ratio of PowerFect™ (µL):DNA (µg) varies from 1:1 to 3:1. We recommend using PowerFect™ (µL):DNA (µg) at 2:1 at a starting point.

The following protocol is given for transfection in 24-well plates, refer to [Table 1](#) for transfection in other culture formats.

- For each well, dilute 0.5 µg of DNA into 50 µl of PowerFect™ Transfection Buffer working solution (1x ) prepared from [Step I](#). Mix by vortexing.
- Add 1.0 µl of PowerFect™ reagent, vortex briefly to mix.
- Incubate for ~10 min at RT to allow PowerFect™/DNA complex to form.

- Note:** Never keep the PowerFect™/DNA complex longer than 20 min.
- Add the PowerFect™/DNA transfection mix to the cells in serum containing medium drop wise.
  - Swirl plate gently to homogenize.
  - Check transfection efficiency 24 to 48 hours post transfection. 48 hours usually give better efficiency.

**Table 1. Recommended Amounts for Different Culture Vessel Formats**

Culture Dish	Culture Medium (mL)	Plasmid DNA (µg)	PowerFect™ Transfection Buffer (1x ) (µL)	PowerFect™ Reagent (µL)
96-well	0.1	0.1	5	0.2
48-well	0.25	0.25	25	0.5
24-well	0.5	0.5	50	1
6-well	2	2.0	200	4
35 mm dish	2	2.0	200	4
60 mm dish	4	4.0	400	8
10 cm / T75	10	10	800	20
15 cm / T175	20	20	1600	40

**Storage:** PowerFect™ siRNA Transfection Reagent is stable for up to 12 months at 4 °C. Always keep PowerFect™ Transfection Buffer (5x ) at RT. This item is shipped at ambient temperature.

## PowerFect™ In Vitro DNA and siRNA Transfection Kit

----- A General Protocol for Transfecting Mammalian Cell

- 100 µl
- 500 µl
- 1000 µl



9601 Medical Center Drive  
Rockville MD 20850  
FAX. 301-560-4919  
TEL. 301-330-5966  
Toll Free. 1-(866)-918-6812  
Email: [info@signagen.com](mailto:info@signagen.com)  
Web: [www.signagen.com](http://www.signagen.com)

This product is for laboratory research ONLY and not for diagnostic use

### Part II. A General Protocol for siRNA Transfection.

#### Step I. Preparation of Working Solution of PowerFect™ Transfection Buffer

PowerFect™ Transfection Buffer (5x ) is provided as 5x concentrated stock solution. To make working solution, dilute one part of the stock solution with 4 parts of ddH<sub>2</sub>O. The 1x PowerFect™ Transfection Buffer is stable at RT for 24 months.

**Note:** Always keep PowerFect™ Transfection Buffer (5x ) at RT. If refrigerated, white precipitates may appear. It won't affect the transfection efficiency. After dilution with 4 parts of ddH<sub>2</sub>O to make PowerFect™ Transfection Buffer (1x ) working solution, the white precipitates will disappear. Always keep PowerFect™ Transfection Buffer working solution (1x ) at RT.

#### Step II. Cell Seeding:

Cells should be plated 18 to 24 hours prior to transfection so that the monolayer cell density reaches to the optimal ~50% confluency at time of transfection. Complete culture medium with serum and antibiotics is freshly added to each well 30-60 min before transfection.

**Note:** High serum levels (>5%) with antibiotics usually do not have inhibitory effect on transfection efficiency. We recommend using complete serum/antibiotics-containing medium to grow the cells during transfection.

#### Step III. Preparation of PowerFect™-siRNA Complex and Transfection Procedures:

For optimal siRNA-mediated silencing, we recommend using 10-80 nM siRNA (final concentration). The following protocol is given for transfection in 6-well plate, refer to [Table 2](#) for transfection in other culture formats.

- For each well, dilute 20 - 160 pmoles siRNA (for a final concentration of 10 to 80 nM per well) into 200 µl of PowerFect™ Transfection Buffer (1x ) prepared from [Step I](#). Mix gently.
- Add 4 µl of PowerFect™ reagent, vortex briefly to mix.
- Incubate for ~10 min at RT to allow PowerFect™/siRNA complexes to form.

- Note:** Never keep the PowerFect™/siRNA complex longer than 20 min.
- Add the PowerFect™/siRNA transfection mix to the cells in serum-containing medium drop wise.
- Swirl plate gently to homogenize.
- Check siRNA silencing efficiency 24 to 72 hours post transfection. 48-72 hours usually give better efficiency.

**Table 2. Recommended Amounts for Different Culture Vessel Formats**

Culture Dish	Culture Medium (mL)	siRNA (pmoles) 10-80 nM	PowerFect™ Transfection Buffer (1x ) (µL)	PowerFect™ Reagent (µL)
96-well	0.1	1 - 8	10	0.4
48-well	0.25	2.5 - 20	25	1
24-well	0.5	5 - 40	50	2
6-well	2	20 - 160	200	4
35 mm dish	2	20 - 160	200	4
60 mm dish	4	40 - 320	400	8
10 cm / T75	10	100 - 800	800	20

**Storage:** PowerFect™ siRNA Transfection Reagent is stable for up to 12 months at 4 °C. Always keep PowerFect™ Transfection Buffer (5x ) at RT. This item is shipped at ambient temperature.